

Junwei Medical

Dental oil free air compressor . High quality air compressor series

- Equipped with improved dry&cooling system, ensuring the supplied air being dry, oil-free, cool and clean. (Optional)
- Special coating process, guaranteeing anti corrosion and antibacterial.
- Copper coil, ensure long working life.
- Set up quickly and convenient.
- Stable and safe.



JW-034A

JW dental air compressor . Dry . Oil free . First choice for cli



JW-031C



JW-032C



JW-036A



JW-032B



JW-032A



JW-038A

Technical Parameters

Unit model Parameters	JW-031C	JW-032B	JW-032A	JW-032C	JW-034A	JW-036A	JW-038A	JW-038A-I
Available units [set]	1	1-2	2-3	3-4	4-5	6-8	9-10	12-15
Voltage [V/Hz]	220/50	220/50	220/50	220/50	220/50	220(380)/50	220(380)/50	380/50
Unit power [W]	550	850	1100	1700	2550	3400	4250	7500
Work pressure [Mpa]	0.5-0.8	0.5-0.8	0.5-0.8	0.5-0.8	0.5-0.8	0.5-0.8	0.5-0.8	0.5-0.8
Flow [L/Min]	62	100	124	200	320	450	560	950
Storage volume [L]	30	45	55	60	90	110	150	180
Noise [db(A)]	≤60	≤62	≤65	≤65	≤68	≤70	≤72	≤75
Product weight [kg]	25	28	43	45	70	92	128	170
Product size [mm]	400*400*510	400*400*620	650*400*600	650*400*620	900*420*670	1154*430*710	1345*450*750	1568*450*825

* The color of real products may have slightly difference from the catalog. Rights reserved for products modifying and designs improving, without further noticing.

