

Junwei-dental

Junwei Medical

Dental suction system . Environment guardians of dental treatment

- Applying of variable frequency technology, domestic initiative technology.
- Two stop separation system, separate air and wast water thoroughly.
- Quickly reach top power, ensuring multiple units work at the same time.
- Stable average suction power on every unit.
- Work safe and continuously.

Essential equipment for high quality dental treatment



Fluctuation pressure sensor:

The suction machine can adjust the power of motor by real-time monitoring the pipeline vacuum degree, thus to guarantee the suction on every unit being used.



Variable frequency control panel:

Real-time monitoring the work situation of the suction machine.



Filter component:

Advanced designed filter component can effectively and practically prevent the inhalation of secretion and blood foam getting in.

Technical Parameters

		JW-X01	JW-Y01	JW-X03	JW-W05
Available units	Set	1	2~3	3~5	5~8
Voltage	V	220~240	220~240	220~240	220~240
Power Frequency	Hz	50/60	50/60	50/60	50/60
Frequency conversion function		X	X	√	√
Unit power	KW	0.4/0.5	0.7/0.8	1.3/1.5	1.5/1.7
Electric current	A	3.0/3.5	4.1/4.8	6.6/6.9	7.6/8.2
Suction power	L/min	≥800	≥1000	≥1100	≥1500
Vacuum	mbar	-110	-130	-180	-180
Noise	dB(A)	≤56	≤56	≤58	≤60
Weight	Kg	18	25	36	38
Size	mm	310*450*585	310*450*585	385*485*615	410*550*615

* The color of real products may have slightly difference from the catalog. Rights reserved for products modifying and designs improving, without further noticing.

Suzhou Junwei Medical equipment Co. Ltd.
 Add: Block 91, No.116 Wutaiشان Road, tech-zone,
 New district, Suzhou
 Tel: 0512 - 67221771 / 67221881
 Fax: 0512 - 67221661
 Website: www.junwei-dental.com
 Official WeChat: Please scan the QR Code on the right.

